



An analytical method for quantifying diazinon, chlorfenapyr, and lufenuron in Napa cabbage using QuEChERS and HPLC-UV

Ahmed Gotah¹, Coralia V. Garcia¹, Sam-Pin Lee^{1,2}, Key Whang^{1*}

¹Department of Food Science and Technology, Keimyung University, Daegu 42601, Korea

²The center for Traditional Microorganism Resource (TMR), Keimyung University, Daegu 42601, Korea

Abstract

A QuEChERS method was developed for the analysis of diazinon, chlorfenapyr, and lufenuron in Napa cabbage. These pesticides represent three different chemical classes and are commonly used in cabbage production in Korea. The objective of the proposed method is a fast, inexpensive, and easy extraction of pesticides, followed by rapid analysis. The proposed method involves a microscale extraction using acetonitrile and dispersive solid phase extraction (SPE), allowing for time and materials savings. The pesticides were separated and quantified using reversed-phase HPLC-UV at 220 nm. The calibration curves showed good linearity ($R^2 > 0.97$), and the limits of detection and quantification were ≤ 0.05 and 1 mg/kg, respectively. Intraday and interday recoveries were in the range 97-116% and 101-112% with $RSD\% \leq 9\%$ for concentrations between 0.5-5 mg/kg. Abnormal recoveries and a substantial matrix effect were initially observed for lufenuron, signaling that optimization of lufenuron recovery requires a slight modification of the method. The proposed method was tested on cabbage samples sold at local markets, which showed no detectable residues of the target pesticides. The proposed method could thus be used for monitoring these pesticides in cabbage and similar vegetables.

Key words : Napa cabbage, pesticide, QuEChERS, HPLC, residue analysis

Introduction

Pesticides are chemicals used in agriculture to combat pests. Although the use of pesticides can provide benefits such as increased inventory and quality of fruits and vegetables, their excessive use has negative environmental effects and can also affect human health (1). Because of this, different national governments have set maximum residue levels (MRLs) for pesticides in various food products.

Napa cabbage, also known as Chinese cabbage, is a common route of exposure of East Asian populations to pesticides. Various insecticides and fungicides are used in Napa cabbage farming, including diazinon, chlorfenapyr, and lufenuron.

Diazinon (O,O-diethyl O-[6-methyl-2-(1-methylethyl)-4-

pyrimidinyl] phosphorothioate) is an organophosphate insecticide-acaricide with contact, stomach, and respiratory action that works by inhibiting cholinesterase (2). The MRL of diazinon as set by the Environmental Protection Agency (EPA) of the United States for brassicas is 0.7 mg/kg (3).

Chlorfenapyr (4-bromo-2-(4-chlorophenyl)-1-ethoxymethyl-5-(trifluoromethyl) pyrrole-3-carbonitrile) is a broad spectrum pyrrole insecticide-acaricide with the appearance of a white-to-tan powder. Chlorfenapyr is actually a pro-insecticide that is converted into an active toxin when ingested by pests, and appears to have little effect on pest predators (4). The MRL of chlorfenapyr as set by the EPA is 1 mg/kg (3).

Lufenuron (N-((2,5-dichloro-4-(1,1,2,3,3,3-hexafluoropropoxy) phenyl) carbamoyl)-2,6-difluorobenzamide) is a benzoylphenylurea insecticide that inhibits chitin synthesis in insects and is used for agricultural and veterinary applications (5). In some countries, lufenuron is only allowed for certain applications; it is allowed in the US in termite baits but not in crops (3).

Multiresidue pesticide analysis is challenging because different pesticides have different solubilities, polarities, and volatilities. Such methods generally require several partitioning

*Corresponding author. E-mail : kwhang@kmu.ac.kr

Phone : 82-53-580-5542, Fax : 82-53-580-5729

Received 10 April 2018; Revised 11 July 2018; Accepted 20 July 2018.

Copyright © The Korean Society of Food Preservation. All rights reserved.

steps followed by chromatography, and as such they are time-consuming and require a substantial volume of solvents. Therefore, the quick, easy, cheap, effective, rugged, and safe (QuEChERS) method was developed in the early 2000s to meet this need (6). The QuEChERS method is a microscale extraction in which pesticides are isolated using a minimal amount of acetonitrile and the extract is cleaned by dispersive solid phase extraction (d-SPE). As such, this method allows for considerable savings of both time and materials, while providing high and reproducible recoveries for a wide range of pesticides in various food products (7). The obtained extract can then be analyzed by HPLC-UV, HPLC-FLV, GC-MS or LC-MS. The advantage of using HPLC is that polar and thermally labile compounds can be analyzed without derivatization, whereas GC analysis may require such a procedure in certain cases. Moreover, HPLC methods tend to be more robust than GC-based ones, and HPLC equipment is more widespread and less expensive than LC-MS (4,8).

The aim of this study was to develop a modified QuEChERS method for quantifying diazinon, chlorfenapyr, and lufenuron residues in Napa cabbage, using HPLC-UV. These three pesticides were selected because they represent three chemical classes of pesticides that are commonly used in Napa cabbage in Korea. HPLC-UV was selected as a low cost and widespread analytical technique suitable for determining most compounds.

Materials and methods

Plant material

Pesticide-free (organic) Napa cabbages were purchased from local shops in Daegu, Korea, and used for validation and recovery experiments. Conventional Napa cabbages were purchased from local markets in Daegu, and analyzed to determine the presence of pesticides.

Chemicals

All solvents used were HPLC-grade. Methanol was obtained from Fisher Scientific (Seoul, Korea) and acetonitrile was obtained from Daejung Chemicals (Siheung, Korea). The water used was from a Milli-Q water purification system (Millipore, Bedford, MA, USA).

Primary secondary amine (PSA) was purchased from Agilent (Little Falls, DE, USA), graphitized carbon black (GCB) and NaCl were obtained from Sigma-Aldrich Co. (St. Louis, MO, USA), and anhydrous MgSO₄ was obtained from

Duksan Pure Chemicals (Ansan, Korea).

Pesticide analytical standards of diazinon (98.5%), chlorfenapyr (purity 99.6%), and lufenuron (purity 99.7%) were purchased from Sigma-Aldrich Co..

Preparation of solutions

Stock solutions containing all analytes (1 mg/mL of each) were prepared in acetonitrile and stored in the dark at -9°C. Serial dilutions in acetonitrile were prepared to make working standard solutions at concentrations of 0.001-1 mg/mL.

QuEChERS samples for calibration

Organic cabbage leaves were chopped, and 200 g of pooled sample was homogenized in a Waring blender (BL233866, Tefal, Shaoxing, China). A QuEChERS method (6) was applied with modifications. Cabbage homogenate (10 g) was weighed in a 50 mL Teflon centrifuge tube, 10 mL of acetonitrile was added, and the tube was shaken for 1 min. Pesticide standards were added to obtain final concentrations of 0.05, 0.1, 0.5, 1.0, 1.5, 2.0, 2.5, and 3.0 mg/kg. Subsequently, the sample tubes were refrigerated for 30 min, 4 g of anhydrous MgSO₄ and 1 g of NaCl were added, and the tubes were shaken for 1 min and centrifuged at 4,000 rpm and 4°C for 5 min. An aliquot of 6 mL of the supernatant was transferred into a 15 mL Teflon centrifuge tube containing 150 mg PSA, 50 mg GCB and 600 mg MgSO₄, and subsequently shaken for 1 min and centrifuged at 4,000 rpm and 4°C for 5 min. An aliquot of 1 mL of the supernatant was then transferred into an Eppendorf tube, and 20 mL of this sample was injected into the HPLC instrument.

Calibration samples were also prepared in solvent in the same way, replacing the blank cabbage with acetonitrile. Reported values represent the averages of triplicate measurements.

QuEChERS samples for recovery calculation

The recovery of the pesticides was determined by spiking cabbage homogenate (10 g) at concentrations of 0.5, 1, and 5 mg/kg. The spiked homogenates were stored in a refrigerator for 3 h to allow pesticide absorption prior to QuEChERS extraction, which was performed as previously described. Five replicates were used to determine recoveries.

Real market samples

Conventional cabbage samples were chosen randomly from markets in the Daegu area and analyzed to determine the presence of the target pesticides. To prepare the samples for

analysis, cabbage leaves were chopped, frozen using liquid nitrogen, and extracted applying the QuEChERS method as previously described.

HPLC analysis

A Shimadzu LC-20A HPLC (Kyoto, Japan) equipped with a UV-Vis detector and a Prominence autosampler was used for analysis. Data were recorded using the LC solution software (Shimadzu). The column was a Waters C₁₈ (3.9×300 mm, 10 μm; Milford, MA, USA).

The mobile phase was methanol:water (75:25, v/v), pumped in isocratic mode at a rate of 1 mL/min for 20 min. The injection volume was 20 μL. The column was maintained at 30°C and detection was performed at 220 nm.

Method validation

The parameters determined were linearity, limit of detection (LOD), limit of quantification (LOQ), precision, accuracy, and recovery, applying previously published equations (9,10).

To determine linearity, the concentrations of the pesticides were determined and the peak areas and heights of the standards were recorded. Calibration equations were obtained by linear regression of the response versus concentration as $y=ax+b$, where y is the response, a is the slope, x is the concentration, and b is the intercept. Calibration curves were obtained for both pesticides diluted in acetonitrile and matrix-matched standards prepared using cabbage blanks for evaluating the matrix effect. The matrix effect was evaluated applying Eq. 1 (11):

$$ME = \frac{\text{slope of matrix-matched calibration curve}}{\text{slope of solvent calibration curve}} \quad (\text{Eq. 1})$$

In this study, the LOD and LOQ were defined as twice the baseline noise and 10× the baseline noise, respectively, in a part of the chromatogram close to the retention time of each analyte (9).

The accuracy of the method was evaluated in terms of recovery, which was calculated using Eq. 2 (10):

$$\%RE = \frac{\text{area of pre-extraction spike}}{\text{area of post-extraction spike}} \times 100 \quad (\text{Eq. 2})$$

The intraday and interday precision of the method was evaluated in terms of repeatability (%RSD), and was calculated using Eq. 3:

$$\%RSD = \frac{\text{STDEV area of post extraction spike}}{\text{AVERAGE area of post extraction spike}} \times 100 \quad (\text{Eq. 3})$$

Results and discussion

QuEChERS and HPLC

The compounds investigated in this study are commonly used pesticides in cabbage farming. Because Koreans frequently eat Napa cabbage (such as in kimchi), cabbage intake may be an important route of exposure to these pesticides in this population. Napa cabbage is also a common component of other East Asian cuisines including those of China and Japan, so this vegetable can also be an important route of pesticide exposure in other Asian populations.

The presence of natural pigments such as chlorophyll complicates the extraction of samples for pesticide analysis. In this study, the use of GCB allowed for the elimination of these interfering pigments, leading to acceptable recoveries. Complete separation of all pesticides was achieved using an isocratic elution with methanol:water (75:25, v/v) as the mobile phase on a C₁₈ column. The separation provided clear resolution after an analysis time of 20 min. Maximum absorption of the pesticides was attained at 220 nm in the UV detector.

The chromatogram of the spiked cabbage sample versus a blank cabbage extract is shown in Fig. 1. The organic cabbage blank provided a clean matrix with no interfering peaks.

Some studies have reported the application of QuEChERS to determine pesticide contents in cabbage, however, these reported methods require the use of mass spectrometers, making them more expensive and requiring more thorough sample preparation. Zhao et al. (12) used LC-MS/MS to determine 183 pesticides in Chinese cabbage, reporting recoveries in the range of 70-120% with LODs of 0.05-3.06 μg/kg. Nguyen et al. (13) used GC-MS to determine 107 pesticides in cabbage, attaining recoveries of 80-115% with LOQs as low as 0.002 mg/kg.

Linearity

Chromatographic data analysis involves peak size measurements, which can be achieved using either peak areas or peak heights. Peak areas are more commonly used; when inadequate resolution could be an issue, like in trace analysis, peak height measurements result in less integration error and would be preferable (14). Thus, for constructing the




Fig. 1. HPLC-UV chromatograms of (A) unspiked organic cabbage and (B) organic cabbage sample spiked with 1 mg/kg of pesticides, registered at 220 nm.

calibration curves, pesticide concentrations of 0.05-3 mg/kg were used and both peak areas and peak heights were integrated. The calibration curves were constructed using standards dissolved in solvent and blank cabbage matrix (Table 1).

Calibration curves obtained using standards dissolved in acetonitrile exhibited higher correlation coefficients (R^2) than those obtained in blank cabbage matrix, particularly when using peak heights for quantification. Calibration curves

obtained in solvent showed good linearity (>0.99) when using peak heights, however, R^2 decreased when using peak areas or blank cabbage matrix (it was higher than 0.97 in all cases). These results suggest interference from the cabbage matrix, requiring further evaluation of the matrix effect.

Matrix interference

The matrix effect can appear as suppression or

enhancement of the analytical signal because of co-eluting matrix components. While this is a well-known problem in LC-MS method development because the ionization process can be easily disrupted by such compounds, it can also occur in conventional techniques such as HPLC-UV and HPLC-FLV. Because of the matrix effect, an analytical signal obtained in solvent may differ from the same signal obtained in a blank matrix (15,16). The matrix effect can be calculated as the ratio between the peak areas corresponding to the post-extraction spike and solvent standard, expressed as a percentage (10) and it can also be evaluated by comparing the slopes obtained in the calibration curves using matrix-matched standards against those obtained with standards diluted in solvent (11).

Limits of detection and quantification

The LOD was defined as twice the signal-to-noise ratio and the LOQ was defined as 10× the signal-to-noise ratio (9). The LODs and LOQs were calculated in matrix-matched standards and are summarized in Table 2. The LODs were 0.03 mg/kg for chlorfenapyr and 0.05 mg/kg for both diazinon and lufenuron. The LOQs were 0.67 mg/kg for chlorfenapyr and 1 mg/kg for both diazinon and lufenuron. For the safety of consumers, the Korean government has set maximum residue limits (MRLs) of 0.05, 0.2, and 0.07 mg/kg for diazinon, chlorfenapyr, and lufenuron respectively in Napa cabbage (18). In Japan, the MRLs for diazinon, chlorfenapyr, and lufenuron in Napa cabbage are 0.1, 2, and 1 mg/kg, respectively (19). Thus, the LODs obtained are below the

Table 1. Calibration data of each pesticide (0.05-3 mg/kg) in acetonitrile and blank cabbage matrix

Compound	Matrix	Height		Area	
		Calibration equation	R ²	Calibration equation	R ²
Diazinon	Acetonitrile	y=1,136.5x+67.857	0.9922	y=22,784x+541.57	0.9892
	Cabbage	y=1,178.5x+119.66	0.9958	y=22,593x+1,760.6	0.9724
Chlorfenapyr	Acetonitrile	y=3,093.4x-59.887	0.9928	y=8,6126x-8,829.3	0.9704
	Cabbage	y=2,943.2x-76.234	0.981	y=73,817x-4,687.9	0.9703
Lufenuron	Acetonitrile	y=2,027.7x-8.9784	0.9927	y=6,5785x-1,944.7	0.9902
	Cabbage	y=999.57x+5.7098	0.9886	y=34,660x+2,661.2	0.9728

Here, the matrix effect was calculated at a concentration of 1 mg/kg using the slopes of the calibration curves obtained by integrating peak heights, and the values were 1.04, 0.95, and 0.49 for diazinon, chlorfenapyr, and lufenuron respectively (Table 2). A value of 1 indicates that there is no matrix effect, whereas values >1 indicate enrichment and values <1 indicate suppression of the analytical signal by the matrix components (11). A matrix effect between 0.9-1.1 can be considered negligible (17). Thus, diazinon and chlorfenapyr exhibited a negligible matrix effect, whereas lufenuron exhibited substantial suppression of the analytical signal.

Therefore, we recommend the use of matrix-matched standards for quantification purposes to counteract the observed matrix effect.

Table 2. LODs, LOQs, and matrix effects of the analyzed pesticides

Compound	LOD (mg/kg)	LOQ (mg/kg)	Matrix effect
Diazinon	0.05	1	1.04
Chlorfenapyr	0.03	0.67	0.95
Lufenuron	0.05	1	0.49

MRLs set in the Korean and Japanese regulations, as well as those set by the EPA for diazinon and chlorfenapyr, although they are higher than those reported by mass spectrometry-based methods (12,13).

Recovery

The recoveries obtained for cabbage samples spiked with

Table 3. Recoveries, intraday and interday precision (% RSD) of the analyzed pesticides

Compound	Concentration (mg/kg)	Recovery 1 st day (%)	Recovery 2 nd day (%)	RSD 1 st day (%)	RSD 2 nd day (%)
Diazinon	5	105	105	8.8	6.8
	1	101	104	8.4	7.7
	0.5	107	104	8.7	8.9
Chlorfenapyr	5	116	112	4.1	1.4
	1	99	101	2.0	4.0
	0.5	103	105	5.0	4.6
Lufenuron	5	199	213	13.7	12.5
	1	129	141	18.0	8.1
	0.5	177	180	13.4	10.6

pesticides at three concentration levels (n=5) are shown in Table 3. Diazinon and chlorfenapyr exhibited intraday recoveries of 101-107% and 99-116% respectively, whereas their interday recoveries were both in the range of 101-112%. The intra and interday %RSDs for these two pesticides were lower than 9%. The acceptable recoveries for trace analysis are 70-120% with %RSD \pm 20% (16), meaning that the values obtained for these two pesticides are adequate. However, in the case of lufenuron, intraday recoveries were in the range of 129-199%, which increased for interday recoveries to 141-213%, even though the %RSD was below 20% in all cases. This result may be due to chromatographic interferences, and suggests that a different clean-up procedure is necessary for the analysis of this particular pesticide (20).

Optimization of lufenuron recovery

Because abnormal recoveries were obtained for lufenuron, the method was modified to improve the recovery of this particular pesticide. Omitting GCB in the clean-up procedure resulted in better recoveries for this pesticide, although more chlorophyll remained in the extract compared to the one processed with GCB. Nevertheless, when analyzing the presence of this particular pesticide, the presence of chlorophyll is of no concern as it does not interfere with the detection of lufenuron.

GCB has been reported to affect the recovery of planar pesticides (21). Thus, omitting GCB may be necessary when analyzing extracts containing planar pesticides or benzoylphenylureas like lufenuron. The recoveries of lufenuron obtained when omitting GCB are shown in Table 4, and are in the range considered optimal; intraday and interday recoveries were 97-116% and 89-106%, respectively.

Table 4. Recoveries, intraday and interday precision (%RSD) of lufenuron in the QuEChERS extract without GCB

Compound	Concentration (mg/kg)	Recovery 1 st day (%)	Recovery 2 nd day (%)	RSD 1 st day (%)	RSD 2 nd day (%)
Lufenuron	5	97	94	1.8	4.2
	1	99	89	4.4	3.1
	0.5	116	106	8.9	9.0

Analysis of real samples

The proposed method was applied to the analysis of 9 cabbage samples selected randomly from markets in Daegu, Korea. No residues of the target pesticides were detected in these samples, indicating that the cabbages analyzed were either free from the presence of these compounds or their

concentrations were below the LOD of the proposed method.

In conclusion, the proposed method allows for simple and rapid determination of diazinon, chlorfenapyr, and lufenuron in Napa cabbage, using QuEChERS and HPLC-UV. This method allows for fast, easy, and inexpensive pesticide analysis without the need for derivatization. This method could be applied to other leaf vegetables for the analysis of various pesticides. However, if planar pesticides or lufenuron are targets, omitting GCB in the extraction process is advisable.

Acknowledgments

This research was supported by the Bisa Research Grant of Keimyung University in 2015.

References

1. Wilson C, Tisdell C (2001) Why farmers continue to use pesticides despite environmental, health and sustainability costs. *Ecol Econ*, 39, 449-462
2. Nemeth-Konda L, Füleky G, Morovjan G, Csokan P (2002) Sorption behaviour of acetochlor, atrazine, carbendazim, diazinon, imidacloprid and isoproturon on Hungarian agricultural soil. *Chemosphere*, 48, 545-552
3. US EPA. About pesticide tolerances. <https://www.epa.gov/pesticide-tolerances/about-pesticide-tolerances> (accessed July 2016)
4. Cao Y, Chen J, Wang Y, Liang J, Chen L, Lu Y (2005) HPLC/UV analysis of chlorfenapyr residues in cabbage and soil to study the dynamics of different formulations. *Sci Total Environ*, 350, 38-46
5. Khay S, Choi JH, Abd El-Aty AM, Mamun MIR, Park BJ, Goudah A, Shin HC, Shim JH (2008) Dissipation behavior of lufenuron, benzoylphenylurea insecticide, in/on Chinese cabbage applied by foliar spraying under greenhouse conditions. *Bull Environ Contam Toxicol*, 81, 369-372
6. Anastassiades M, Lehotay S, Stajnbaher D, Schenk FJ (2003) Fast and easy multiresidue method employing acetonitrile extraction/partitioning and "dispersive solid-phase extraction" for the determination of pesticide residues in produce. *J AOAC Int*, 86, 412-431
7. Wilkowska A, Biziuk M (2011) Determination of pesticide residues in food matrices using the QuEChERS

- methodology. *Food Chem*, 125, 803-812
8. Soodi M, Garshasbi A, Eskandari S (2012) Development a HPLC method for simultaneous determination of Azinphose methyl, Diazinon, Phosalone and Chlorpyrifos residues in fruit. *J Pharm Health Sci*, 1, 79-87
 9. Francotte E, Davatz A, Richert P (1996) Development and validation of chiral high-performance liquid chromatography methods for the quantitation of valsartan and of the tosylate of valinebenzyl ester. *J Chromatogr B*, 686, 77-83
 10. Kruve A, Kunnapas A, Herodes K, Leito I (2008) Matrix effects in pesticide multi-residue analysis by liquid chromatography-mass spectrometry. *J Chromatogr A*, 1187, 58-66
 11. Romero-Gonzalez R, Garrido Frenich A, Martinez Vidal JL, Prestes OD, Grió SL (2011) Simultaneous determination of pesticides, biopesticides and mycotoxins in organic products applying a quick, easy, cheap, effective, rugged and safe extraction procedure and ultra-high performance liquid chromatography-tandem mass spectrometry. *J Chromatogr A*, 1218, 1477-1485
 12. Zhao MA, Feng YN, Zhu YZ, Kim JH (2014) Multi-residue method for determination of 238 pesticides in Chinese cabbage and cucumber by liquid chromatography-tandem mass spectrometry: Comparison of different purification procedures. *J Agric Food Chem*, 62, 11449-11456
 13. Nguyen TD, Yu JE, Lee DM, Lee GH (2008) A multiresidue method for the determination of 107 pesticides in cabbage and radish using QuEChERS sample preparation method and gas chromatography mass spectrometry. *Food Chem*, 110, 207-213
 14. Niec D, Kunicki PK (2015) Validation of an assay for quantification of free normetanephrine, metanephrine and methoxytyramine in plasma by high performance liquid chromatography with coulometric detection: Comparison of peak-area vs. peak-height measurements. *J Chromatogr B*, 1002, 63-70
 15. Matuszewski BK, Constanzer ML, Chavez-Eng CM (2003) Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC-MS/MS. *Anal Chem*, 75, 3019-3030
 16. Sampaio MRF, Tomasini D, Cardoso LV, Caldas SS, Primel EG (2012) Determination of pesticide residues in sugarcane honey by QuEChERS and liquid chromatography. *J Braz Chem Soc*, 23, 197-205
 17. Dong B, Zhao Q, Hu J (2015) Dissipation kinetics of emamectin benzoate and lufenuron residues in cabbage grown under field conditions. *Environ Monit Assess*, 187, 765
 18. Korean Ministry of Food and Drug Safety. Pesticides and veterinary drugs information. <http://www.foodnara.go.kr/residue/main.do> (accessed July 2016)
 19. The Japan Food Chemical Research Foundation. Maximum residue limits (MRLs) list of agricultural chemicals in foods. <http://www.m5.ws001.squarestart.ne.jp/foundation/search.html> (accessed July 2016)
 20. Li M, Liu X, Dong F, Xu J, Kong Z, Li Y, Zheng Y (2013) Simultaneous determination of cyflumetofen and its main metabolite residues in samples of plant and animal origin using multi-walled carbon nanotubes in dispersive solid-phase extraction and ultrahigh performance liquid-chromatography-tandem mass spectrometry. *J Chromatogr A*, 1300, 95-103
 21. Akiyama Y, Matsuoka T, Mitsuhashi T (2009) Multi-residue screening method of acidic pesticides in agricultural products by liquid chromatography/time of flight mass spectrometry. *J Pestic Sci*, 34, 265-272